

# Comparison of test results for several analytes using capillary blood versus venous blood on the Atellica CH and Atellica IM Analyzers

J. Freeman,<sup>1</sup> J. Cheek,<sup>1</sup> S. Singh,<sup>1</sup> C. DiPasquale,<sup>2</sup> E.G. Olson,<sup>2</sup> R. Barr,<sup>2</sup> H. McKenzie,<sup>2</sup> J.W. Jacobson<sup>2</sup>

<sup>1</sup>Siemens Healthcare Diagnostics Inc., Tarrytown, NY, USA; <sup>2</sup>Babson Diagnostics, Inc., Austin, TX, USA

## Background

Capillary blood collection is a minimally invasive method used inside and outside of traditional healthcare settings, especially important for home monitoring, pharmacies, and those with difficult veins, or needle phobia. Capillary blood is suitable for diagnostic testing at the point of care and within core laboratories. This study compared results from capillary\* and venous blood samples for common chemistry (CH) and immunoassay (IM) methods on core laboratory Atellica CH and Atellica IM Analyzers.<sup>1,21</sup>

## Methods

- An Institutional Review Board-approved clinical study was performed at Babson Diagnostics, Inc. with ambulatory adult subjects recruited with informed consent.
- The specimen equivalence study was designed per CLSI EP35.<sup>22</sup> Matched set samples were collected: venous blood using Vacutainer SST and EDTA (comparators) and capillary blood using Microtainer SST and MAP EDTA tubes. A total of 97–115 matched sets were assayed. Samples were collected on a rolling basis and assayed in duplicate within 4 hours of collection using one reagent lot on one Atellica CH or Atellica IM Analyzer.
- Data were analyzed by regression of the capillary vs venous comparator samples, utilizing the first replicate from each sample for slope, intercept, and correlation coefficient. Relative repeatability was estimated by pooled SD ratio with 95% CI for each parameter.
- Chemistry assays included Enzymatic Hemoglobin A1c (A1c\_E), Alanine Aminotransferase (ALT), Albumin (AlbP), Alkaline Phosphatase (ALP\_2c), Aspartate Aminotransferase (AST), Calcium (Arsenazo III) (CA\_2), Calcium (cresolphthalein) (Ca), Total Cholesterol (Chol\_2), electrolytes (A-LYTE IMT Na K Cl), Diazo Direct Bilirubin (D\_DBil), Diazo Total Bilirubin (D\_TBil), Enzymatic Creatinine (ECre3), Glucose Hexokinase (GluH\_3), HDL Cholesterol (HDLc), LDL Cholesterol (LDLc), Total Protein (TP\_2), Triglycerides (Trig\_2), Urea nitrogen (UN\_c).
- Immunoassays included Free T4 (FT4), Total T3 (T3), and Thyroid Stimulating Hormone (TSH3-ULII).

## Results


Assay slopes ranged from 0.92–1.04 and correlation coefficients (r) were >0.950, except for Potassium (0.936). Overall, testing covered at least 80% of the entire assay range. The repeatability SD ratio for capillary to venous samples ranged from 0.59–1.33, with ratio 95% CI lower bounds ranging from 0.41–0.92 and upper bounds ranging from 0.86–1.91.

The design requirements for specimen equivalence were met for capillary versus venous samples on the Atellica CH and IM Analyzers. When analyzed by Deming regression, Weighted Deming regression, or Passing Bablok the capillary samples on the Atellica CH and IM Analyzers recovered samples spanning at least 80% of the measuring interval, with slopes of  $1.00 \pm 0.10$  and a correlation coefficient (r)  $\geq 0.950$  (except for Potassium) compared to venous samples on the Atellica CH and IM Analyzers. Figure 1 presents Deming, Passing Bablok fit, and difference plots for select assays listed in Table 1, for capillary vs venous samples on the Atellica CH or IM Analyzers.


Table 1. Specimen equivalency results for capillary (Y) versus venous (X) serum

| Analyte                           | Assay               | Slope | Intercept    | Sample Interval     | R value | n   | SD, Y/X (95% CI) |
|-----------------------------------|---------------------|-------|--------------|---------------------|---------|-----|------------------|
| Whole Blood                       |                     |       |              |                     |         |     |                  |
| Hemoglobin A1c                    | A1c_E               | 1.00  | +0.03 %NGSP  | 4.5–13.32 %NGSP     | 0.998   | 109 | 0.79 (0.54,1.15) |
| Serum                             |                     |       |              |                     |         |     |                  |
| Alanine Aminotransferase          | ALT                 | 0.99  | +0 U/L       | 7–871 U/L           | 1.000   | 114 | 1.00 (0.69,1.45) |
| Albumin                           | AlbP                | 1.00  | +0.0 g/dL    | 3.1–7.5 g/dL        | 0.992   | 114 | 1.00 (0.69,1.45) |
| Alkaline Phosphatase              | ALP_2c              | 1.00  | +0 U/L       | 40–883 U/L          | 1.000   | 114 | 0.78 (0.49,1.03) |
| Aspartate Aminotransferase        | AST                 | 0.99  | +2 U/L       | 11–847 U/L          | 1.000   | 109 | 0.78 (0.53,1.13) |
| Calcium                           | CA_2 (Arsenazo III) | 1.00  | +0.0 mg/dL   | 4.5–13.4 mg/dL      | 0.983   | 114 | 0.86 (0.59,1.24) |
| Calcium                           | Ca                  | 1.00  | +0.0 mg/dL   | 4.2–14.1 mg/dL      | 0.983   | 114 | 0.63 (0.43,0.9)  |
| Total Cholesterol                 | Chol_2              | 1.00  | +0 mg/dL     | 74–581 mg/dL        | 0.999   | 115 | 1.00 (0.69,1.44) |
| Chloride                          | A-LYTE IMT Cl       | 1.00  | +2.00 mEq/L  | 78.0–166 mEq/L      | 0.990   | 114 | 1.09 (0.75,1.58) |
| Potassium                         | A-LYTE IMT K        | 0.92  | +0.49 mEq/L  | 2.01–8.98 mEq/L     | 0.936   | 99  | 1.00 (0.67,1.49) |
| Sodium                            | A-LYTE IMT Na       | 1.00  | +0 mEq/L     | 70.0–197 mEq/L      | 0.991   | 115 | 1.04 (0.72,1.51) |
| Direct Bilirubin                  | D_DBil (Diazo)      | 0.93  | -0.01 mg/dL  | 0.10–10.71 mg/dL    | 1.000   | 97  | 0.91 (0.61,1.36) |
| Total Bilirubin                   | D_TBil (Diazo)      | 0.98  | -0.01 mg/dL  | 0.14–19.86 mg/dL    | 1.000   | 115 | 0.67 (0.46,0.96) |
| Creatinine                        | ECre3 (Enzymatic)   | 0.98  | -0.03 mg/dL  | 0.34–25.88 mg/dL    | 1.000   | 114 | 0.90 (0.62,1.30) |
| Glucose                           | GluH_3 (Hexokinase) | 1.03  | +3 mg/dL     | 48–649 mg/dL        | 0.995   | 114 | 1.00 (0.69,1.45) |
| HDL Cholesterol                   | HDLc                | 1.00  | -0.0 mg/dL   | 29.0–178.9 mg/dL    | 0.999   | 115 | 1.33 (0.92,1.91) |
| LDL Cholesterol                   | LDLc                | 0.99  | +1 mg/dL     | 32–359 mg/dL        | 0.998   | 115 | 0.84 (0.58,1.21) |
| Total Protein                     | TP_2                | 1.03  | -0.2 g/dL    | 3.9–10.6 g/dL       | 0.990   | 115 | 1.00 (0.69,1.44) |
| Triglycerides                     | Trig_2              | 1.00  | +3 mg/dL     | 27–904 mg/dL        | 0.999   | 115 | 0.91 (0.63,1.31) |
| Urea Nitrogen                     | UN_c                | 1.00  | +0 mg/dL     | 7–134 mg/dL         | 0.999   | 115 | 1.20 (0.83,1.73) |
| Free Thyroxine (Free T4)          | FT4                 | 1.04  | +0.0 ng/dL   | 0.8–10.6 ng/dL      | 0.999   | 115 | 1.00 (0.69,1.44) |
| Total Triiodothyronine (Total T3) | T3                  | 0.98  | +0.02 ng/mL  | 0.50–7.20 ng/mL     | 0.993   | 115 | 1.02 (0.70,1.47) |
| Thyroid Stimulating Hormone (TSH) | TSH3-ULII           | 1.00  | -0.002 mIU/L | 0.128–135.638 mIU/L | 1.000   | 111 | 0.59 (0.41,0.86) |


A. Atellica CH Enzymatic Hemoglobin A1c (A1c\_E) assay




B. Atellica CH Alkaline Phosphatase (ALP\_2c) assay




C. Atellica CH Total Cholesterol (Chol\_2) assay




D. Atellica CH Diazo Total Bilirubin (D\_TBil) assay




E. Atellica CH Enzymatic Creatinine (ECre3) assay




F. Atellica CH Total Protein (TP\_2) assay




G. Atellica CH Urea Nitrogen (UN\_c) assay



H. Atellica IM Free Thyroxine (FT4) assay



I. Atellica IM Total Triiodothyronine (T3) assay



J. Atellica IM Thyroid Stimulating Hormone (TSH3-ULII) assay




Figure 1. Representative Deming and Passing Bablok regression and difference plots for capillary serum compared to venous serum specimen equivalence on the Atellica CH and IM analyzers for several of the assays in Table 1.

## Conclusion

Results for routinely tested Atellica CH and Atellica IM assays were equivalent between capillary and venous samples across the assay range measured.

## References

- Atellica CH Alanine Aminotransferase (ALT) assay, 11110168\_EN Rev. 06, 2024-09.
- Atellica CH Albumin (AlbP) assay, 11642849\_EN Rev. 03, 2022-05.
- Atellica CH Alkaline Phosphatase, Concentrated (ALP) assay, 11110163\_EN Rev. 03, 2019-07.
- Atellica CH Aspartate Aminotransferase (AST) assay, 11642857\_EN Rev. 04, 2024-10.
- Atellica CH Calcium, 2 (Arsenazo III) (CA\_2) assay, 11116511\_EN Rev. 03, 2019-03.
- Atellica CH Calcium (Ca) cresoph assay, 11110158\_EN Rev. 02, 2019-10.
- Atellica CH Cholesterol, 2 (Chol\_2) assay, 11110172\_EN Rev. 04, 2020-11.
- Atellica Diazo Direct Bilirubin (D\_DBil) (D\_DBil) assay, 11537291\_EN Rev. 02, 2023-09.
- Atellica CH Diazo Total Bilirubin (D\_TBil) assay, 11537295\_EN Rev. 02, 2023-11.
- LYTE Integrated Multisensor (IMT Na K Cl) assay, 11109447\_EN Rev. 06, 2023-10.
- Atellica CH Enzymatic Creatinine 3 (ECre3) assay, 11537279\_EN Rev. 04, 2024-11.
- Atellica CH Glucose Hexokinase, 3 (GluH\_3) assay, 11110155\_EN Rev. 06, 2025-04.
- Atellica CH HDL Cholesterol (HDLc) assay, 11537569\_EN Rev. 05, 2024-08.
- Atellica CH LDL Cholesterol (LDLc) assay, 11537570\_EN Rev. 06, 2024-08.
- Atellica CH Triglycerides, 2 (Trig\_2) assay, 11537302\_EN Rev. 04, 2022-03.
- Atellica CH Urea Nitrogen (UN\_c) assay, 11110156\_EN Rev. 03, 2020-08.
- Atellica CH Total Protein (TP\_2) assay, 11110167\_EN Rev. 03, 2020-08.
- Atellica IM Free T4 (FT4) assay, 10995348\_EN Rev. 07, 2023-03.
- Atellica IM T3 Total (T3) assay, 10995424\_EN Rev. 04, 2023-03.
- Atellica IM Thyroid Stimulating Hormone (TSH3-ULII) assay, 11202198\_EN Rev. 05, 2023-03.
- Atellica CH Enzymatic Hemoglobin A1c (A1c\_E) assay, 11645641\_EN Rev. 04, 2025-03.
- CLSI. Assessment of Equivalence or Suitability of Specimen Types for Medical Laboratory Measurement Procedures. 1st ed. CLSI guideline EP35. Wayne, PA: Clinical and Laboratory Standards Institute; 2019.

\*Under development. Not available for sale in the U.S.A.

All trademarks are the property of their respective owners

For educational and scientific exchange only.

Not for sales or promotional use.

Published by Siemens Healthcare Diagnostics Inc.

© Siemens Healthcare Diagnostics Inc., 2026

Scan QR code for  
downloadable copy of this  
poster or visit  
www.medicalaffairs.  
siemens-healthineers.com

