

M. Borgman, D. Blagovcanin, Z. Eischeid, A. Sinopoli, J. Cheek, N. Ozgen, J. Snyder
Siemens Healthcare Diagnostics Inc., Tarrytown, NY, U.S.

Total bilirubin (TBil), a byproduct of hemoglobin catabolism, is a critical biomarker in the evaluation of neonatal hyperbilirubinemia. In newborns, even modest elevations in bilirubin levels are closely monitored due to the risk of kernicterus, a serious and potentially irreversible neurological complication associated with excessive bilirubin exposure. Accurate measurement of TBil is therefore essential for early identification, risk stratification, and clinical management of hyperbilirubinemia in neonates. This study evaluated the reference range and method comparison of the Diazo TBil (D₅ TBil) assay¹ on the Atellica CH analyzer using prospectively collected neonatal serum samples.

- Healthy neonatal serum samples with no history of hyperbilirubinemia were prospectively collected (23 samples: 0–1 day; 27 samples: 1–2 days; 26 samples: 3–5 days).
- All samples were tested using the Atellica CH D₂Tbil assay.
- Results were compared with published reference ranges, and ranges were considered verified if $\leq 10\%$ of the samples fell outside of the established limits.

- A total of 113 prospectively collected serum samples from healthy neonates and those diagnosed with hyperbilirubinemia, along with 11 contrived samples prepared by spiking neonatal serum with unconjugated bilirubin, were tested using the Atellica CH 930 D_TBil and Roche COBAS 8000 TBIL (BILT3) assays.
- Weighted Deming regression was used to assess slope, y-intercept, and correlation between the two methods.
- Method comparison studies were conducted consistent with the governing standard *CLSI EP09C-ED3: Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. Individual neonate serum samples were tested on the Atellica CH system and Roche COBAS 8000 system. The method comparison was completed with the Atellica CH D_TBil assay as the test assay and the Roche COBAS 8000 TBIL assay as the predicate assay. Individual samples were obtained from outside vendors through the external testing site and specimen acquisition group. Eleven samples were spiked with unconjugated bilirubin material supplied with Certificate of Analysis from Frontier Specialty Chemicals to cover the assay measuring interval. A single replicate was processed for each sample on both systems.

- The Atellica CH D_Tbil assay demonstrated complete concordance with physician-assigned classifications of healthy versus hyperbilirubinemic neonates when evaluated against the established clinical nomogram.²
- Age specific expected Tbil values were verified for healthy neonates aged 0–1 day (<8.00 mg/dL), 1–2 days (<12.00 mg/dL), and 3–5 days (<16.00 mg/dL), consistent with clinical expectation.

Specimen Type	Neonate Age (Day)	Reference Interval (SI)	# of Samples	# of Samples > RI ^a	% Samples > RI
Serum	0–1	< 8.00 mg/dL (< 136.80 µmol/L)	23	0	0
Serum	1–2	< 12.00 mg/dL (< 205.20 µmol/L)	27	0	0
Serum	3–5	< 16.00 mg/dL (< 273.60 µmol/L)	26	0	0

The verification testing data for healthy neonates were acceptable and met the acceptance criteria, consistent with clinical expectation.

Neonate Age (Day)	Reference Interval
0–1	< 8.00 mg/dL (< 136.80 μmol/L)
1–2	< 12.00 mg/dL (< 205.20 μmol/L)
3–5	< 16.00 mg/dL (< 273.60 μmol/L)

The Atellica CH D_TBil assay demonstrated agreement with the Roche COBAS TBIL assay across the neonatal concentration range, 0.10–21.81 mg/dL, with a regression of $y = 1.06x - 0.02$ mg/dL and $r = 0.997$, indicating strong analytical concordance within clinically relevant bilirubin levels.

Sample Type	Comparison Assay	Comparison Analyzer (x)	n	r	Regression Equation	Sample Range
Serum	Roche COBAS TBIL (BILT3)	Roche COBAS 8000	124	0.997	$y = 1.06x - 0.02 \text{ mg/dL}$ ($y = 1.06x - 0.34 \text{ }\mu\text{mol/L}$)	0.10–21.81 (1.71–372.95)

The design requirements for method comparison were met for the Atellica CH D_Tbil assay compared to the Roche COBAS TBIL (BILT3) assay. When analyzed by Weighted Deming regression, the two measurement methods recovered samples spanning the measuring interval, with a slope of 1.00 ± 0.10 , intercept 0.00 ± 0.10 mg/dL, and a correlation coefficient (r) ≥ 0.950 .

Left Plot: Comparison of Atellica CH 930 D_TBil Assay and Roche COBAS 8000 TBIL Assay

Y-axis: Atellica CH 930 D_TBil Assay (mg/dL)
X-axis: Roche COBAS 8000 TBIL Assay (mg/dL)

Regression equation: $y = 1.06x - 0.02$ mg/dL
Correlation coefficient: $r = 0.997$
Sample size: $n = 124$

Legend: Δ Data, — Weighted Deming fit

Right Plot: Difference between Atellica CH 930 D_TBil Assay and Roche COBAS 8000 TBIL Assay

Y-axis: Difference (%)
X-axis: Roche COBAS 8000 TBIL Assay (mg/dL)

Using neonatal serum samples, the Atellica CH D_T Bil assay demonstrated clinically acceptable agreement with the Roche COBAS TBIL (BILT3) assay and verified age-appropriate neonatal reference values, supporting its use for accurate assessment and management of neonatal hyperbilirubinemia.

Published by Siemens Healthcare Diagnostics Inc.
QR700049117
© Siemens Healthcare Diagnostics Inc., 2026



Scan QR code for
downloadable copy of this
poster or visit
[www.medicalaffairs.
siemens-healthineers.com](http://www.medicalaffairs.siemens-healthineers.com)