

Development of pre-analytics for a Rapid, Integrated Sepsis Diagnostic and Prognostic Workflow combining Pathogen Detection with Host Immune Response Profiling

Divya Sharma Khandige¹, Josh Chenoweth², Ramya VM¹, Deeksha Hebbar¹, Thomas A Hall³, Lisa Risen³, Christine Marzan³, Michael Mosel³, James Hannis³, Hila Roshanravan⁴, and Megan Carpenter⁴

¹Siemens Healthcare Private Limited, Bangalore, IN. ²The Henry M. Jackson Foundation for the Advancement of Military Medicine, Inc., Bethesda, Maryland, USA;

³Janus I-Science, San Diego, USA; ⁴Siemens Healthcare Diagnostics, Tarrytown, USA.

Introduction

Sepsis is a life-threatening syndrome caused by a dysregulated immune response to infection and remains a major global health burden, contributing to approximately 20% of deaths worldwide¹. Early identification of both the infecting pathogen and illness severity is critical, particularly in austere and military medicine settings where timely treatment or transport decisions are essential. However, the current gold standard, blood culture, requires 0.5 to 5 days for pathogen identification and antimicrobial susceptibility testing, delaying clinical decision-making^{2,3,4}. No commercially available integrated solution currently combines pathogen detection with sepsis severity assessment. To address this gap, we developed INSPIRE (INtegrated Sepsis Pathogen identification and Immune Response Evaluation), a 6-hour nanopore sequencing–based diagnostic and prognostic workflow for direct whole-blood analysis (Figure 1). INSPIRE enables unbiased pathogen identification, antimicrobial resistance (AMR) detection, and mortality risk prediction using universally amplified DNA signatures and relative gene expression profiling of 18 clinically relevant host mRNAs⁵, including three reference genes.




Figure 1. Schematic representation of the overview of INSPIRE workflow.

Methods

Nucleic acid purification and processing

Host mRNA was enriched from PAXgene blood samples using an optimized capture probe-based method targeting 18 transcripts. The 18-mRNA panel was identified in an earlier study and described by Chenoweth et al. 2024. Pathogen DNA was recovered from EDTA whole blood using a universal lysis method compatible with bacterial and fungal pathogens, followed by capture probe-based purification. DNA and RNA eluates were amplified by multiplex PCR. Primer sets for pathogen identification were designed to target short, highly conserved genomic regions flanking systematic phylogenetic variation across a broad range of bacterial and fungal species. In case of host mRNA conserved exon-exon junction was used to design PCR primers. Amplicons were purified using AMPure beads, prepared for nanopore sequencing and run on ONT flow cells.

Software and algorithms

The real-time analysis of this integrated sequencing run was done using “NanoSepsIDx” software, and reads were mapped against the proprietary database to identify pathogen/s and its AMR. In addition, a sepsis severity risk score for the same patient sample is determined by the integrated “Predict_d28” classifier algorithm that uses a weighted fuzzy-logic approach to predict mortality based on sepsis endotypes.

Integrated workflow for pathogen identification and risk prediction

Whole blood samples collected in EDTA tubes were spiked with Staphylococcus aureus or Enterococcus faecium at a final concentration of 50 CFU/mL. These contrived EDTA samples were processed in parallel with matched PAXgene whole blood samples for host mRNA purification. Pathogen DNA was purified using the capture probe-based method, whereas host mRNA was purified using either the capture probe-based method or the standard silica-based method as a control. A total of 10 sample sets were processed to evaluate the integrated workflow using the capture probe–based method for both pathogen identification and risk prediction, as well as the combined approach of capture probe-based pathogen identification and silica-based host mRNA purification for risk prediction. In this preliminary workflow format, the data is continuously captured, analyzed and reports are generated at 5 different time points of 15, 30, 60, 120, and 180 minutes.

Results

Sequencing summary results from four healthy donors showed that the capture probe–based mRNA purification method yielded results comparable to those obtained with the standard silica-based total RNA purification method, followed by sequencing. Samples purified using the capture probe–based method were sequenced on the same Oxford Nanopore Technologies (ONT) flow cell as total RNA. Each healthy donor sample was processed in duplicates using the capture probe assay, along with one standard sample, and the results are summarized in Table 1.

Table 1. Sequencing summary results for four healthy donors comparing the capture probe–based mRNA purification method with the standard silica-based total RNA purification method used for template preparation.

	Sample 1			Sample 2			Sample 3			Sample 4		
Sample	CP Rep 1	CP Rep 2	Standard	CP Rep 1	CP Rep 2	Standard	CP Rep 1	CP Rep 2	Standard	CP Rep 1	CP Rep 2	Standard
28-day mortality risk assessment:	Low	Low	Low	Low	Low	Low	Low	Low	Low	–	–	Low
28-day mortality risk weight:	0.069	0.098	0.069	0.046	0.06	0.116	0.053	0.083	0.142	–	–	0.122

Note: CP: capture probe-based mRNA purification. Standard: Total RNA purified using silica columns.

Interday and intraday repeatability of the capture probe-based mRNA purification method was evaluated using four healthy donor samples processed in duplicate across three independent days and analyzed by nanopore sequencing. Replicates from each donor were subsequently analyzed within a single sequencing run. Overall, the results demonstrated consistent performance across both intraday and interday conditions, supporting the repeatability of the capture probe–based mRNA purification workflow.




Figure 2. mRNA read counts normalized to the reference mRNA read count, plotted on a logarithmic scale on the y-axis. The x-axis represents the 15 different mRNAs analyzed in the study. The four panels present results from four different healthy donors. The similar expression patterns observed across donors indicate consistent mRNA profiling and support the repeatability of the capture probe–based purification and sequencing workflow.

Sequencing generated a compiled report that included pathogen ID, AMR gene detection, and sepsis risk prediction. AMR detection results are not shown, as the strains used to spike the samples did not contain resistance genes. Detection of AMR genes was evaluated and confirmed in a separate study.

